


Product Datasheet

WIBG (PYM1) antibody (orb1330177)

Description	Mouse monoclonal antibody to WIBG (PYM1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human WIBG (NP_115721) produced in HEK293T cell.
Target	PYM1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q9BRP8
NCBI	NM_032345
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...