

Product Datasheet

Peripherin (PRPH) antibody (orb1330275)

Description Mouse monoclonal antibody to Peripherin (PRPH)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human PRPH (NP_006253) produced in HEK293T cell.

Target PRPH

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 53.5 kDa

Uniprot ID [P41219](#)

NCBI [NM_006262](#)

Dilution Range WB 1:2000, IHC 1:150

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...