


Product Datasheet

DNAAF11 antibody (orb1330296)

Description	Mouse monoclonal antibody to DNAAF11
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human LRRC6 (NP_036604) produced in HEK293T cell.
Target	DNAAF11
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	54.1 kDa
Uniprot ID	Q86X45
NCBI	NM_012472
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...