


Product Datasheet

SUOX antibody (orb1330398)

Description	Mouse monoclonal antibody to SUOX
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 320-545 of human SUOX(NP_000447) produced in E.coli.
Target	SUOX
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	60.1 kDa
Uniprot ID	P51687
NCBI	NM_000456
Dilution Range	WB 1:500, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...