


Product Datasheet

Rb (RB1) antibody (orb1330470)

Description	Mouse monoclonal antibody to Rb (RB1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 309-590 of human RB1(NP_000312) produced in E.coli.
Target	RB1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P06400
NCBI	NM_000321
Dilution Range	WB 1:200, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...