

Product Datasheet

Estrogen Related Receptor beta (ESRRB) antibody (orb1330479)

Description Mouse monoclonal antibody to Estrogen Related Receptor beta (ESRRB)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Human recombinant protein fragment corresponding to amino acids 1-252 of human ESRRB (NP_004443) produced in E.coli.

Target ESRRB

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2b

Clonality Monoclonal


MW 56 kDa

Uniprot ID [O95718](#)

NCBI [NM_004452](#)

Dilution Range WB 1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...