

Product Datasheet

TFPI2 antibody (orb1330544)

Description	Mouse monoclonal antibody to TFPI2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 23-235 of human TFPI2 (NP_006519) produced in E.coli.
Target	TFPI2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	26.8 kDa
Uniprot ID	P48307
NCBI	NM_006528
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

