


Product Datasheet

p16INK4A (CDKN2A) antibody (orb1330582)

Description	Mouse monoclonal antibody to p16INK4A (CDKN2A)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Full length human recombinant protein of human CDKN2A(NP_000068) produced in E.coli.
Target	CDKN2A
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	16.4 kDa
Uniprot ID	Q8N726
NCBI	NM_000077
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...