

Product Datasheet

Constitutive androstane receptor (NR1I3) antibody (orb1330597)

Description	Mouse monoclonal antibody to Constitutive androstane receptor (NR1I3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-236 of human NR1I3 (NP_001070950) produced in E.coli.
Target	NR1I3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	40.3 kDa
Uniprot ID	Q14994
NCBI	NM_001077482
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

