


Product Datasheet

uPA (PLAU) antibody (orb1330621)

Description	Mouse monoclonal antibody to uPA (PLAU)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 107-379 of human PLAU (NP_002649) produced in E.coli.
Target	PLAU
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	46.3 kDa
Uniprot ID	P00749
NCBI	NM_002658
Dilution Range	WB 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...