

Product Datasheet

p16INK4A (CDKN2A) antibody (orb1330667)

Description Mouse monoclonal antibody to p16INK4A (CDKN2A)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen Full length human recombinant protein of human CDKN2A (NP_000068) produced in E.coli.

Target CDKN2A

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 16.4 kDa

Uniprot ID [Q8N726](#)

NCBI [NM_000077](#)

Dilution Range WB 1:2000, IHC 1:150, IF 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...