


Product Datasheet

Progesterone Receptor (PGR) antibody (orb1330687)

Description	Mouse monoclonal antibody to Progesterone Receptor (PGR)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-298 of human PGR(NP_000917) produced in E.coli.
Target	PGR
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	98.8 kDa
Uniprot ID	P06401
NCBI	NM_000926
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...