


Product Datasheet

Nkx3.1 (NKX3-1) antibody (orb1330695)

Description	Mouse monoclonal antibody to Nkx3.1 (NKX3-1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human NKX3-1 (NP_006158) produced in E.coli.
Target	NKX3-1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	26.2 kDa
Uniprot ID	Q99801
NCBI	NM_006167
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...