


Product Datasheet

CD229 antibody (orb13307)

Description	CD229 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CD229 (301-400/655 aa)
Target	CD229
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67 kDa
Uniprot ID	Q9HBG7
Expiration Date	12 months from date of receipt.



Paraffin-
embedded
human lung
carcinoma f...