

Product Datasheet

Apolipoprotein E (APOE) antibody (orb1330724)

Description Mouse monoclonal antibody to Apolipoprotein E (APOE)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC

Immunogen Full length human recombinant protein of human APOE (NP_000032) produced in HEK293T cell.

Target APOE

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 34.2 kDa

Uniprot ID [P02649](#)

NCBI [NM_000041](#)

Dilution Range IHC 1:150

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...