

Product Datasheet

Thyroid Hormone Receptor alpha (THRA) antibody (orb1330733)

Description Mouse monoclonal antibody to Thyroid Hormone Receptor alpha (THRA)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Human recombinant protein fragment corresponding to amino acids 1-220 of human THRA (NP_955366) produced in E.coli.

Target THRA

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2a


Clonality Monoclonal

Uniprot ID [P10827](#)

NCBI [NM_199334](#)

Dilution Range WB 1:2000

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...