


Product Datasheet

Lipin 3 (LPIN3) antibody (orb1330770)

Description	Mouse monoclonal antibody to Lipin 3 (LPIN3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human LPIN3 (NP_075047) produced in HEK293T cell.
Target	LPIN3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	93.4 kDa
Uniprot ID	Q9BQK8
NCBI	NM_022896
Dilution Range	WB 1:500, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...