

Product Datasheet

Galectin 3 (LGALS3) antibody (orb1330909)

Description Mouse monoclonal antibody to Galectin 3 (LGALS3)

Species/Host Mouse

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Full length human recombinant protein of human LGALS3 (NP_002297) produced in E.coli.

Target LGALS3

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG2a

Clonality Monoclonal


MW 26 kDa

Uniprot ID [P17931](#)

NCBI [NM_002306](#)

Dilution Range WB 1:2000, IHC 1:150

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...