

Product Datasheet

SGK196 (POMK) antibody (orb1330934)

Description	Mouse monoclonal antibody to SGK196 (POMK)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 141-350 of human SGK196 (NP_115613) produced in E.coli.
Target	POMK
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	39.9 kDa
Uniprot ID	Q9H5K3
NCBI	NM_032237
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.

