

Product Datasheet

Laminin 5 (LAMB3) antibody (orb1331019)

Description	Mouse monoclonal antibody to Laminin 5 (LAMB3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 424-773 of human LAMB3 (NP_000219) produced in E.coli.
Target	LAMB3
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	127.6 kDa
Uniprot ID	Q13751
NCBI	NM_000228
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

