


Product Datasheet

S100 alpha 6 (S100A6) antibody (orb1331076)

Description	Mouse monoclonal antibody to S100 alpha 6 (S100A6)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human S100A6 (NP_055439) produced in E.coli.
Target	S100A6
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	10 kDa
Uniprot ID	P06703
NCBI	NM_014624
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...