


Product Datasheet

gamma IL2 Receptor antibody (orb13311)

Description	gamma IL2 Receptor antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	KLH conjugated synthetic peptide derived from human IL-2R gamma (51-150/369 aa)
Target	IL-2R gamma
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	P34902
Dilution Range	ICC=1:100, Flow-Cyt=1µg /test
Expiration Date	12 months from date of receipt.



ICC analysis of SHSY5Y cell using IL-2R ...