


Product Datasheet

N4BP2L2 antibody (orb1331133)

Description	Mouse monoclonal antibody to N4BP2L2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human N4BP2L2 (NP_055702) produced in E.coli.
Target	N4BP2L2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	67.3 kDa
Uniprot ID	Q92802
NCBI	NM_014887
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...