


Product Datasheet

PPP2R5D antibody (orb1331194)

Description	Mouse monoclonal antibody to PPP2R5D
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human PPP2R5D (NP_006236) produced in E.coli.
Target	PPP2R5D
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	69.8 kDa
Uniprot ID	Q14738
NCBI	NM_006245
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...