


Product Datasheet

SETD2 antibody (orb1331205)

Description	Mouse monoclonal antibody to SETD2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1787-2144 of human SETD2 (NP_054878) produced in E.coli.
Target	SETD2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9BYW2
NCBI	NM_014159
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of Human recombinant...