



Product Datasheet

SLAMF9 antibody (orb13313)


Description	SLAMF9 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human SLAMF9 (51-150/289 aa)
Target	SLAMF9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30 kDa
Uniprot ID	Q9D780
Dilution Range	WB=1:500-2000, Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Blank control:
(mo)splenocytes(blue)
Iso...



Sample: A549
(human)cell Lysate at
40 ug...



Blank control (black
line): Raji (black)...