

Product Datasheet

CK1 epsilon (CSNK1E) antibody (orb1331332)

Description	Mouse monoclonal antibody to CK1 epsilon (CSNK1E)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 153-416 of human CSNK1E (NP_689407) produced in E.coli.
Target	CSNK1E
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	47.1 kDa
Uniprot ID	P49674
NCBI	NM_152221
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.

