


Product Datasheet

Von Willebrand Factor (VWF) antibody (orb1331473)

Description	Mouse monoclonal antibody to Von Willebrand Factor (VWF)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 764-1263 of human VWF (NP_000543) produced in SF9 cell.
Target	VWF
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	309.3 kDa
Uniprot ID	P04275
NCBI	NM_000552
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...