


Product Datasheet

IL1RA (IL1RN) antibody (orb1331482)

Description	Mouse monoclonal antibody to IL1RA (IL1RN)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human IL1RN (NP_776215) produced in E.coli.
Target	IL1RN
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	16 kDa
Uniprot ID	P18510
NCBI	NM_173843
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...