


Product Datasheet

TTF1 (NKX2-1) antibody (orb1331598)

Description	Mouse monoclonal antibody to TTF1 (NKX2-1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full length human recombinant protein of human NKX2-1 (NP_003308) produced in E.coli.
Target	NKX2-1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	38.4 kDa
Uniprot ID	P43699
NCBI	NM_003317
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...