

Product Datasheet

CD35 antibody (orb13316)

Description

CD35 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human CD35. Please contact us for the exact immunogen sequence. The peptide is available as orb374760.

Target

CD35. Recognized endogenous levels of total CD35 antibody.

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

224 kDa

Uniprot ID

P17927

NCBI

21, 000573, 31, 000564

Entrez


1378

Dilution Range


WB: 1:100-2000, IHC-P: 1:100-800 (based on 0.5 mg/ml)

Expiration Date


12 months from date of receipt.



Western blot analysis of rat lymph node ...



Immunohistochemical staining of paraffin...



IHC-P image of mouse lymph node tissue u...