

Product Datasheet

Nogo B receptor (NUS1) antibody (orb1331663)

Description	Mouse monoclonal antibody to Nogo B receptor (NUS1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 38-116 of human NUS1 (NP_612468) produced in E.coli.
Target	NUS1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q96E22
NCBI	NM_138459
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

