


Product Datasheet

MRPS11 antibody (orb1331794)

Description	Mouse monoclonal antibody to MRPS11
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human MRPS11 (NP_073750) produced in E.coli.
Target	MRPS11
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	20.4 kDa
Uniprot ID	P82912
NCBI	NM_022839
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...