

Product Datasheet

Macrophage Scavenger Receptor I (MSR1) antibody (orb1331847)

Description	Mouse monoclonal antibody to Macrophage Scavenger Receptor I (MSR1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 197-451 of human MSR1(NP_619730) produced in E.coli.
Target	MSR1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	42.8 kDa
Uniprot ID	P21757
NCBI	NM_138716
Dilution Range	WB 1:200 - 1:1000
Expiration Date	12 months from date of receipt.

