


Product Datasheet

CD44v6 antibody (orb13319)

Description	CD44v6 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from rat CD44v6 (325-429/429 aa)
Target	CD44
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53 kDa
Uniprot ID	P16070
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...