


Product Datasheet

Ki67 (MKI67) antibody (orb1331962)

Description	Mouse monoclonal antibody to Ki67 (MKI67)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1160-1493 of human MKI67 (NP_002408) produced in E.coli.
Target	MKI67
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P46013
NCBI	NM_002417
Dilution Range	IHC 1:100~500, WB 1:200 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (10ug)...