

Product Datasheet

GAPDH antibody (orb1331971)

Description	Mouse monoclonal antibody to GAPDH
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human GAPDH (NP_002037) produced in HEK293T cell.
Target	GAPDH
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	35.9 kDa
Uniprot ID	P04406
NCBI	NM_002046
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.




Figure from citation:
Western blot analy...




Figure from citation:
Western blot analy...




Figure from citation:
Western blot analy...