


Product Datasheet

PI 3 Kinase Class 2A (PIK3C2A) antibody (orb1332310)

Description	Mouse monoclonal antibody to PI 3 Kinase Class 2A (PIK3C2A)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 230-560 of human PIK3C2A (NP_002636) produced in E.coli.
Target	PIK3C2A
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	190.5 kDa
Uniprot ID	O00443
NCBI	NM_002645
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of E.coli lysate (...)