


Product Datasheet

NFH (NEFH) antibody (orb1332539)

Description	Mouse monoclonal antibody to NFH (NEFH)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 778-1020 of human NEFH (NP_066554) produced in E.coli.
Target	NEFH
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P12036
NCBI	NM_021076
Dilution Range	WB 1:2000, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...