


Product Datasheet

Vimentin (VIM) antibody (orb1332552)

Description	Mouse monoclonal antibody to Vimentin (VIM)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 210-466 of human VIM (NP_003371) produced in E.coli.
Target	VIM
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	53.5 kDa
Uniprot ID	P08670
NCBI	NM_003380
Dilution Range	WB 1:200 - 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...