

Product Datasheet

CD69 antibody (orb13326)

Description

CD69 antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human CD69. Please contact us for the exact immunogen sequence. The peptide is available as orb12905.

Target

CD69

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

23 kDa

Uniprot ID

[Q07108](#), [P37217](#)

NCBI

[11](#), [21](#), [001772](#), [001781](#)

Entrez


[969](#)

Dilution Range


IF/ICC: 1:50-400, IHC-P: 1:50-400, WB: 1:200-2000

Expiration Date


12 months from date of receipt.



Western blot analysis of mouse spleen (l...



IHC-P staining of human thyroid cancer t...



Immunohistochemical staining of paraffin...