


Product Datasheet

CTNND1 antibody (orb1332685)

Description	Mouse monoclonal antibody to CTNND1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human CTNND1 (NP_001078927) produced in HEK293T cell.
Target	CTNND1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	108 kDa
Uniprot ID	O60716
NCBI	NM_001085458
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...