


Product Datasheet

Ki67 (MKI67) antibody (orb1332690)

Description	Mouse monoclonal antibody to Ki67 (MKI67)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2924-3256 of human MKI67 (NP_002408) produced in E.coli
Target	MKI67
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P46013
NCBI	NM_002417
Dilution Range	WB 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of E.coli lysate (...)