


Product Datasheet

DDX59 antibody (orb1332842)

Description	Mouse monoclonal antibody to DDX59
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-303 of human DDX59 (NP_001026895) produced in E.coli.
Target	DDX59
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	68.6 kDa
Uniprot ID	Q5T1V6
NCBI	NM_001031725
Dilution Range	IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...