


Product Datasheet

CD4 antibody (orb1332868)

Description	Mouse monoclonal antibody to CD4
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human CD4 (NP_000607) produced in HEK293T cell.
Target	CD4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	48.3 kDa
Uniprot ID	P01730
NCBI	NM_000616
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...