


Product Datasheet

ZFAND5 antibody (orb1332887)

Description	Mouse monoclonal antibody to ZFAND5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZFAND5 (NP_005998) produced in E.coli.
Target	ZFAND5
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	23 kDa
Uniprot ID	O76080
NCBI	NM_006007
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...