

Product Datasheet

CD83 antibody (orb13329)

Description

CD83 antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human CD83. Please contact us for the exact immunogen sequence. The peptide is available as orb12906.

Target

CD83

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW


23 kDa

Uniprot ID
[Q01151](#), [O88324](#)
NCBI
[001035370](#), [001040280](#), [11](#)
Entrez
[9308](#)
Dilution Range


IHC-P: 1:50-300, WB: 1:200-2000

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of rat lymph node ...



Immunohistochemical staining of paraffin...