


Product Datasheet

Isocitrate dehydrogenase (IDH1) antibody (orb1332936)

Description	Mouse monoclonal antibody to Isocitrate dehydrogenase (IDH1)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide around the R132 region of the human IDH conjugated to KLH
Target	IDH1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	46.5 kDa
Uniprot ID	075874
NCBI	NM_005896
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of Equivalent amou...