


Product Datasheet

ARL2BP antibody (orb1332959)

Description	Mouse monoclonal antibody to ARL2BP
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human ARL2BP (NP_036238) produced in E.coli.
Target	ARL2BP
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	18.6 kDa
Uniprot ID	Q9Y2Y0
NCBI	NM_012106
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...