

Product Datasheet

Meis homeobox 3 (MEIS3) antibody (orb1333027)

Description Mouse monoclonal antibody to Meis homeobox 3 (MEIS3)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, WB

Immunogen Human recombinant protein fragment corresponding to amino acids 1-261 of human MEIS3 (NP_064545) produced in E.coli.

Target MEIS3

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 46 kDa

Uniprot ID [Q99687](#)

NCBI [NM_020160](#)

Dilution Range WB 1:2000, IF 1:100, FLOW 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...